

Recombinant Murine TNF-related apoptosis-inducing Ligand/TNFSF10



Catalog No: #AP60179

Orders: order@signalwayantibody.com

Package Size: #AP60179-1 10ug #AP60179-2 100ug #AP60179-3 500ug

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Murine TNF-related apoptosis-inducing Ligand/TNFSF10
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Uniprot	P50592
GeneID	22035
Calculated MW	Approximately 20.2 kDa, a single non-glycosylated polypeptide chain containing 175 amino acids.
Target Sequence	MPRGGRPQKV AAHITGITRR SNSALIPISK DGKTLGQKIE SWESSRKGHS FLNHVLFNRNG ELVIEQEGLY YIYSQTYFRF QEAEDASKMV SKDKVRTKQL VQYIYKYTSY PDPIVLMKSA RNSCWSRDAE YGLYSIQGG LFELKKNDRIFVSVTNEHLM DLDQEASFFG AFLIN
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, with 3 mM DTT.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only