

Recombinant Murine soluble Tumor Necrosis Factor Receptor Type I/TNFRSF1A



Catalog No: #AP60182

Orders: order@signalwayantibody.com

Package Size: #AP60182-1 5ug #AP60182-2 100ug #AP60182-3 500ug

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Murine soluble Tumor Necrosis Factor Receptor Type I/TNFRSF1A
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	P25118
GeneID	21937
Calculated MW	Approximately 21.1 kDa, a single non-glycosylated polypeptide chain containing 191 amino acids.
Target Sequence	IHPSGVTGLV PSLGDREKRD SLCPQGKYVH SKNNSICCTK CHKGTYLVSD CPSPGRDVC RECEKGTFTA SQNYLRQCLS CKTCRKEMSQ VEISPCQADK DTVCGCKENQ FQRYLSETHF QCVDCSPCFN GTVTIPCKET QNTVCNCHAG FFLRESECV CSHCKKNEEC MKLCLPPPLA NVTNPQDSGT A
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only