

Recombinant Murine Nephroblastoma-overexpressed Gene Protein Homolog

Catalog No: #AP60194

Orders: order@signalwayantibody.com

Package Size: #AP60194-1 5ug #AP60194-2 100ug #AP60194-3 500ug

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Murine Nephroblastoma-overexpressed Gene Protein Homolog
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Uniprot	Q64299
GeneID	18133
Calculated MW	Approximately 36.4 kDa, a single non-glycosylated polypeptide chain containing 333 amino acids.
Target Sequence	QVSASLRCPSP RCPPKCPSPS PTCAPGVRSV LDGCSCCPVC ARQRGESCSE MRPCDQSSGL YCDRSADPNN QTGICMVPEG DNCVFDGVIY RNGEKFEFPC QYFCTCRDGG IGCLPRCQLD VLLPGPDCPA PRKVAVPGEC CEKWTGCSDE QGTQGTLLGL ALPAYRPEAT VGVEVSDSSI NCIEQTTEWS ACSKSCGMGV STRVTNRNRQ CEMVKQTRLC IVRPCEQEPE EVTDKKGKKC LRTKSLKAI HLQFENCTSL YTYKPRFCGV CSDGRCCTPH NTKTIQVEFQ CLPGEIHKP VMVIGTCTCY SNCPQNNEAF LQDLELKTSR GEI
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris-HCl, pH 8.6, 150 mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only