

Recombinant Murine Interferon-lambda2/Interleukin-28A

Catalog No: #AP60200


Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP60200-1 5ug #AP60200-2 100ug #AP60200-3 500ug

Description

Product Name	Recombinant Murine Interferon-lambda2/Interleukin-28A
Host Species	E.coli
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Uniprot	Q4VK74
GeneID	330496
Calculated MW	Approximately 19.7 kDa, a single non-glycosylated polypeptide chain containing 174 amino acids.
Target Sequence	DPVPRATRLP VEAKDCHIAQ FKSLSPKELQ AFKAKDAIE KRLLEKDMRC SSSLISRAWD LKQLQVQERP KALQAEVALT LKVWENMTDS ALATILGQPL HTLSHIHSQL QTCTQLQATA EPKPPSRRLS RWLHRLQEAQ SKETPGCLEL SVTSNLFRL TRDLKCVASG DQCV
Formulation	LyophilizedB fromB aB 0.2B umB filteredB concentratedB solutionB inB PBS.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only