

# Recombinant Murine Mesencephalic Astrocyte-Derived Neurotrophic Factor

Catalog No: #AP60206


Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP60206-1 5ug #AP60206-2 100ug #AP60206-3 500ug

## Description

Product Name	Recombinant Murine Mesencephalic Astrocyte-Derived Neurotrophic Factor
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Uniprot	Q9CXI5
Calculated MW	Approximately 18.2 kDa, a single non-glycosylated polypeptide chain containing 158 amino acids.
Target Sequence	LRPGDCEVCI SYLGRFYQDL KDRDVTFSPA TIEEELIKFC REARGKENRL CYYIGATDDA ATKIINEVSK PLAHHIPVEK ICEKLVKKDS QICELKYDNQ IDLSTVDLKK LRVKELKKIL DDWGENCKGC AEKSDYIRKI NELMPKYAPK AASARTDL
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

## Images



Note: This product is for in vitro research use only