

Recombinant Murine Persephin

Catalog No: #AP60210



Package Size: #AP60210-1 5ug #AP60210-2 100ug #AP60210-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------|--|
| Product Name | Recombinant Murine Persephin |
| Host Species | Escherichia coli. |
| Purification | > 95 % by SDS-PAGE and HPLC analyses. |
| Uniprot | O70300 |
| Calculated MW | Approximately 20.7 kDa, a disulfide-linked homodimer of two 96 amino acid polypeptide chains. |
| Target Sequence | ALAGSCLWS LTLPVAELGL GYASEEKVIF RYCAGSCPQE ARTQHSLVLA RLRGRGRAHG RPCCQPTSYA DVTFLDDQHH WQLPQLSAA ACGCGG |
| Formulation | Lyophilized from a 0.2 µm filtered concentrated solution in 30 % ACN, 0.1 % TFA, 150 mM NaCl. |
| Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution. |

Note: This product is for in vitro research use only