

Recombinant Rat Interleukin-9

Catalog No: #AP60233


Package Size: #AP60233-1 10ug #AP60233-2 100ug #AP60233-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Interleukin-9
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	D4A819
GeneID	116558
Calculated MW	Approximately 14.3 kDa, a single non-glycosylated polypeptide chain containing 127 amino acids.
Target Sequence	MQRCSTSWG I QHTSYLIENL KDPSSK CSC SANVT SCLCL PIPSD DCTTP CFQEGMSQVT NATQQSKFSP FFFRVKRIVE TLKSNKCQFF SCEKPCNQTT AGNTVSFLKS LLKTFQKTEV QVQRSRA
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only