

Recombinant Rat Interleukin-10

Catalog No: #AP60234


Package Size: #AP60234-1 10ug #AP60234-2 100ug #AP60234-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Interleukin-10
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	P29456
GeneID	25325
Calculated MW	Approximately 18.6 kDa, a single non-glycosylated polypeptide chain containing 160 amino acids.
Target Sequence	SKGHSIRGDN NCTHFPVSQT HMLRELRAAF SQVKTFQKK DQLDNILLTD SLLQDFKGYL GCQALSEMIK FYLVEVMPQA ENHGPEIKEH LNSLGEKLT LWIQLRRCHR FLPCENKSKA VEQVKNDFNK LQDKGVYKAM NEFDIFINCI EAYVTLKMKN
Formulation	Lyophilized from a 0.2 µm filtered solution in 20 mM Tris-HCl, pH 8.0, 100 mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only