

Recombinant Rat Interleukin-17A

Catalog No: #AP60237



Package Size: #AP60237-1 5ug #AP60237-2 100ug #AP60237-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Interleukin-17A
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Uniprot	Q61453
Calculated MW	Approximately 30.0 kDa, a disulfide-linked homodimer of two 133 amino acid polypeptide chains.
Target Sequence	AVLIPQSSVC PNAEANNFLQ NVKVNKVIN SLSSKASSRR PSDYLNRRSTS PWTLNRNEDP DRYPSVIWEA QCRHQRCVNA EGKLDHHMNS VLIQQEILVL KREPEKCPFT FRVEKMLVGV GCTCVSSIVR HAS
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 2 x PBS, pH 6.5.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only