

Recombinant Rat Leukemia Inhibitory Factor

Catalog No: #AP60245


Package Size: #AP60245-1 5ug #AP60245-2 100ug #AP60245-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Leukemia Inhibitory Factor
Host Species	Escherichia coli.
Purification	> 96 % by SDS-PAGE and HPLC analyses.
Uniprot	P17777
Calculated MW	Approximately 19.8 kDa, a single non-glycosylated polypeptide chain containing 180 amino acids.
Target Sequence	SPLPITPVNA TCAIRHPCHG NLMNQIKSQL AQLNGSANAL FISYYTAQGE PFPNNVDKLC APNMTDFPPF HANGTEKTKL VELYRMVTYL GASLTNITWD QKNLNPTAVS LQIKLNATTD VMRGLLSSVL CRLCNKYHVG HVDVPCVPDN SSKEAFQRKK LGCQLLGTYK QVISVLAQAF
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 1 x PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only