

Recombinant Rat Vascular Endothelial Growth Factor 164, Yeast



Catalog No: #AP60253

Orders: order@signalwayantibody.com

Package Size: #AP60253-1 10ug #AP60253-2 100ug #AP60253-3 500ug

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Vascular Endothelial Growth Factor 164, Yeast
Host Species	Yeast
Purification	> 95 % by SDS-PAGE and 90% by SEC-HPLC analyses.
Calculated MW	Theoretically as a disulfide-linked homodimeric protein, the product consists of two 165 amino acid polypeptide chains. As a result of glycosylation, it migrates to at least two bands with molecular weights ranging from 25.7 kDa in SDS-PAGE under reducing conditions.
Target Sequence	MAPTTEGEQK AHEVVKFMDV YQRSYCRPIE TLVDIFQEYP DEIEYIFKPS CVPLMRCAGC CNDEALECVP TSESNVTMQI MRIKPHQSQH IGEMSFLQHS RCECRPKKDR TKPEKHCEPC SERRKHLFVQ DPQTCKCSCK NTDSRCKARQ LELNERTCRC DKPRR
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only