

## Recombinant Rat Platelet-derived Growth Factor-BB

Catalog No: #AP60255


Package Size: #AP60255-1 10ug #AP60255-2 100ug #AP60255-3 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                 |  |
|-----------------|--|
| Product Name    | Recombinant Rat Platelet-derived Growth Factor-BB  |
| Host Species    | Escherichia coli.  |
| Purification    | > 98 % by SDS-PAGE and HPLC analyses.  |
| Uniprot         | Q05028   |
| GeneID          | 24628  |
| Calculated MW   | Approximately 24.4 kDa, a disulfide-linked homodimeric protein containing two 109 amino acid residues polypeptide (B chain).   |
| Target Sequence | SLGSLAAAEP AVIAECKTRT EVFQISRNLI DRTNANFLVW PPCVEVQRCS GCCNNRNVQC RASQVQMRPV<br>QVRKIEIVRK KPVFKKATVT LEDHLACKCE TVVTPRPVT   |
| Formulation     | Lyophilized from a 0.2 µm filtered concentrated solution in 1 x PBS, pH 6.9.   |
| Storage         | Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution. |

## Images



Note: This product is for in vitro research use only