

Recombinant Rat Glia Maturation Factor beta

Catalog No: #AP60257


Package Size: #AP60257-1 10ug #AP60257-2 100ug #AP60257-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Glia Maturation Factor beta
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Other Names	GMFB
Uniprot	Q63228
GeneID	81661
Calculated MW	Approximately 16.6 kDa, a single non-glycosylated polypeptide chain containing 141 amino acids.
Target Sequence	SESLVVCDVA EDLVEKLRKF RFRKETHNAA IIMKIDKDKR LVLDEELEG VSPDELKDEL PERQPRFIVY SYKYQHDDGR VSYPLCFIFS SPLGCKPEQQ MMYAGSKNKL VQTAELTKVF EIRNTEDLTE EWLREKLGFF H
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Background

The glia maturation factor beta belongs to the actin-binding proteins ADF family, GMF subfamily. It contains an ADF-H domain, but the research of crystallography and NMR reveals that there are structures different between human and mouse ADF-H domain. GMF-β is involved in the differentiation, maintenance, and regeneration of the nervous system. It also inhibition of proliferation of tumor cells.

Note: This product is for in vitro research use only