

Recombinant Rat Resistin

Catalog No: #AP60266



Package Size: #AP60266-1 5ug #AP60266-2 100ug #AP60266-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Resistin
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	Q8K4J7
GeneID	246250
Calculated MW	Approximately 20.2 kDa, a disulfide-linked homodimer consisting of two 95 amino acid polypeptide chains.
Target Sequence	MPSMSLCPMD EAISKINQD FSSLLPAAMK NTVLHCWSVS SRGLASCPE GTTVTSCSCG SGCGSWDVRE DTMCHCQCGS IDWTAARCCT LRVGS
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH7.4, with 0.02 % Tween-20.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only