

Recombinant Porcine Interleukin-1 Receptor Antagonist Protein

Catalog No: #AP60270


Package Size: #AP60270-1 5ug #AP60270-2 100ug #AP60270-3 500ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Porcine Interleukin-1 Receptor Antagonist Protein
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	Q29056
GeneID	397499
Calculated MW	Approximately 17.1 kDa, a single non-glycosylated polypeptide chain containing 152 amino acids.
Target Sequence	HPLGKRPCRMAAFRIWDVNVKTFYLRNNQLVAGYLQGPNTELEEKIDVVPVEPHFVFLGIHGGKLCCLSCVKSGDEMKLQLDAVNITDLRKNSQDKRFTFIRSDSGPTTSFESAACPGWFLCTALEADQPVGLTNTPKAAVKVTKFYFQQDQ
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 1 x PBS, pH 7.4, 1 mM DTT.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only