

## Recombinant Equine Interferon-gamma

Catalog No: #AP60276


Package Size: #AP60276-1 100ug #AP60276-2 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Equine Interferon-gamma
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	P42160
GeneID	100034181
Calculated MW	Approximately 16.7 kDa, a single non-glycosylated polypeptide chain containing 143 amino acids.
Target Sequence	QA AFFKEIEN LKEYFNASNP DVG DGGPLFL DILKNWKEDS DKKIIQSQIV SFYFKLFENL KDNQVIQKSM DTIKEDLFVK FFNSSTSKLE DFQKLIQIPV NDLKVQRKAI SELIKVMNDL SPKANLRKRK RSQNPFRGRR ALQ
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 2 x PBS, pH 7.4, with 5 % trehalose.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

## Images



Note: This product is for in vitro research use only