

# Recombinant Human Granulocyte Chemotactic Protein 2/CXCL6



Catalog No: #AP60284

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP60284-1 5ug #AP60284-2 100ug #AP60284-3 500ug

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Granulocyte Chemotactic Protein 2/CXCL6
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Uniprot	P80162
GeneID	6372
Calculated MW	Approximately 7.9 kDa, a single non-glycosylated polypeptide chain containing 72 amino acids.
Target Sequence	VLTELRCTCL RVTLRVNPKT IGKLQVFPAG PQCSKVEVVA SLKNGKQVCL DPEAPFLKKV IQKILDGSGNK KN
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only