

Recombinant Human Lymphocyte Activation Gene 1 Protein/CCL4L1

Catalog No: #AP60300

Orders: order@signalwayantibody.com


Package Size: #AP60300-1 5ug #AP60300-2 100ug #AP60300-3 500ug

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Lymphocyte Activation Gene 1 Protein/CCL4L1
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	Q8NHW4
GeneID	388372
Calculated MW	Approximately 7.8 kDa, a single non-glycosylated polypeptide chain containing 69 amino acids.
Target Sequence	APMGSDPPTA CCFSYTARKL PRNFVVDYVE TSSLCSQPAV VFQTKRGKQV CADPSESWWQ EYVYDLELN
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only