

Recombinant Rhesus Macaque BCA-1/CXCL13

Catalog No: #AP60324



Package Size: #AP60324-1 5ug #AP60324-2 100ug #AP60324-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rhesus Macaque BCA-1/CXCL13
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Calculated MW	Approximately 10.3 kDa, a single non-glycosylated polypeptide chain containing 87 amino acids.
Target Sequence	VLEVYYTHLR CRCVQESSVF IPRRFIDRIQ ISPRGNGCPR KEIIVWKKNK SVVCVDPQAE WIQRIMEMLR KKSSSTPPVP VFKRKIP
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM Tris-HCl, pH 8.0, 300mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only