

Recombinant Murine LPS-induced CXC Chemokine, 70a.a./CXCL5

Catalog No: #AP60334

Package Size: #AP60334-1 5ug #AP60334-2 100ug #AP60334-3 500ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Murine LPS-induced CXC Chemokine, 70a.a./CXCL5
Host Species	Escherichia coli.
Purification	> 96 % by SDS-PAGE and HPLC analyses.
Uniprot	P50228
GeneID	20311
Calculated MW	Approximately 7.6 kDa, a single non-glycosylated polypeptide chain containing 70 amino acids.
Target Sequence	TELRCVCLTV TPKINPKLIA NLEVIPAGPQ CPTVEVIAKL KNQKEVCLDP EAPVIKKIIQ KILGSDKKKA
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only