

## Recombinant Murine I-TAC/CXCL11

Catalog No: #AP60337



Package Size: #AP60337-1 5ug #AP60337-2 100ug #AP60337-3 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Murine I-TAC/CXCL11
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Uniprot	Q9JHH5
GeneID	56066
Calculated MW	Approximately 9.1 kDa, a single non-glycosylated polypeptide chain containing 79 amino acids.
Target Sequence	FLMFKQGRCL CIGPGMKAVK MAEIEKASVI YPSNGCDKVE VIVTMKAHKR QRCLDPRSKQ ARLIMQAIEK KNFLRRQNM
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 10 mM Sodium Citrate, pH 4.0, with 600 mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only