

Recombinant Murine Lymphotactin/XCL1

Catalog No: #AP60343



Package Size: #AP60343-1 5ug #AP60343-2 100ug #AP60343-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Murine Lymphotactin/XCL1
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	P47993
GeneID	16963
Calculated MW	Approximately 10.1 kDa, a single non-glycosylated polypeptide chain containing 93 amino acids.
Target Sequence	VGTEVLEESS CVNLQTQRLP VQKIKTYIIW EGAMRAVIFV TKRGLKICAD PEAKWVKAAL KTVDGRASTR KNMAETVPTG AQRSTSTAIT LTG
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only