

Recombinant Murine Fractalkine/CX3CL1

Catalog No: #AP60344



Package Size: #AP60344-1 5ug #AP60344-2 100ug #AP60344-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Murine Fractalkine/CX3CL1
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	O35188
GeneID	20312
Calculated MW	Approximately 8.7 kDa, a single non-glycosylated polypeptide chain containing 76 amino acids and comprises only the chemokine domain of Murine Fractalkine.
Target Sequence	QHLGMTKCEI MCGKMTSRIP VALLIRYQLN QESCGKRAIV LETTQHRRFC ADPKEKWWQD AMKHLDHQAA ALTKNG
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only