

# Recombinant Rat Breast and Kidney-expressed Chemokine/CXCL14

Catalog No: #AP60372

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP60372-1 5ug #AP60372-2 100ug #AP60372-3 500ug

## Description

Product Name	Recombinant Rat Breast and Kidney-expressed Chemokine/CXCL14
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Uniprot	Q8K453
GeneID	306748
Calculated MW	Approximately 9.4 kDa, a single non-glycosylated polypeptide chain containing 77 amino acids.
Target Sequence	SKCKCSRKGP KIRYSDVKKL EMKPKYPHCE EKMVIITKS MSRYRQGEHC LHPKLQSTKR FIKWYNAWNE KRRVYEE
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 30 % Acetonitrile and 0.1 % TFA.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only