

Recombinant Rat Macrophage Inflammatory Protein-1 alpha/CCL3



Catalog No: #AP60376

Orders: order@signalwayantibody.com

Package Size: #AP60376-1 5ug #AP60376-2 100ug #AP60376-3 500ug

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Macrophage Inflammatory Protein-1 alpha/CCL3
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Uniprot	P50229
GeneID	25542
Calculated MW	Approximately 7.9 kDa, a single non-glycosylated polypeptide chain containing 69 amino acids.
Target Sequence	APYGADTPTA CCFSYGRQIP RKFIADYFET SSLCSQPGVI FLTKRNRQIC ADPKETWVQE YITELELNA
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 30 % Acetonitrile and 0.1 % TFA.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only