

Recombinant Rat Monocyte Chemoattractant Protein 3/CCL7

Catalog No: #AP60379


Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP60379-1 10ug #AP60379-2 100ug #AP60379-3 500ug

Description

Product Name	Recombinant Rat Monocyte Chemoattractant Protein 3/CCL7
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Other Names	Protein FIC, Small-inducible Cytokine A7, Intercrine/chemokine MARC
Uniprot	Q9QXY8
GeneID	287561
Calculated MW	Approximately 8.5 kDa, a single non-glycosylated polypeptide chain containing 74 amino acids.
Target Sequence	QPDGTSSTC CYVKKQKIPK RNLKSYRKIT SSRCPWEAVI FKTKKGMEVC AEAHQKWVEE AIAYLDMKTS TPKP
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in 2 x PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Background

Rat CCL7, belonging to the CC chemokine family, is also named as MCP-3, because the sequence comparison suggested that the rat CCL7 is a homologue of the human MCP-3 and is postulated the rat MCP-3. The MCP-3 protein family signals through CCR2 expect MCP-1 and possess cross-reacts across species. CCL7/MCP3 has chemotactic functions for monocytes and eosinophils, but not for neutrophils. Additionally, it can activated NK cells as well as CD4+ and CD8+ T lymphocytes. CCL7 was found to interact with MMP2.

Note: This product is for in vitro research use only