

## Recombinant Rat Prolactin

Catalog No: #AP60406


Package Size: #AP60406-1 10ug #AP60406-2 100ug #AP60406-3 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Rat Prolactin
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Uniprot	P01237
GeneID	24683
Calculated MW	Approximately 22.6 kDa, a single non-glycosylated polypeptide chain containing 197 amino acids.
Target Sequence	LPVCSGGDCQ TPLPELFDRV VMLSHYIHTL YTDMFIEFDK QYVQDREFIA KAINDCPTSS LATPEDKEQA QKVPPEVLLN LILSLVHSWN DPLFQLITGL GGIHEAPDAI ISRAKEIEEQ NKRLLEGIEK IISQAYPEAK GNEIYLWSQ LPSLQGVDEE SKDLAFYNNI RCLRRDSHKV DNYLKFLRCQ IVHKNNC
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

## Images



Note: This product is for in vitro research use only