

Recombinant PreScission Protease, Low Endotoxin

Catalog No: #AP60429



Package Size: #AP60429-1 100IU #AP60429-2 250IU

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant PreScission Protease, Low Endotoxin
Host Species	Escherichia coli
Other Names	3C protease, Picornain 3C, PSP
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 6 months from date of receipt, -20 to -70 °C as supplied.- 3 months, -20 to -70 °C under sterile conditions after opening.

Note: This product is for in vitro research use only