

# Recombinant Human Ubiquitin-conjugating Enzyme E2 M

Catalog No: #AP60439

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP60439-1 5ug #AP60439-2 100ug #AP60439-3 500ug

## Description

Product Name	Recombinant Human Ubiquitin-conjugating Enzyme E2 M
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Uniprot	P61081
GeneID	9040
Calculated MW	Approximately 20.9 kDa, a single non-glycosylated polypeptide chain containing 183 amino acids.
Target Sequence	MIKLFSLKQQ KKEEESAGGT KGSSKASAA QLRIQKDINE LNLPKTCDIS FSDPDDLNF KLVICPDEGF YKSGKFVFSF KVGQGYPHDP PKVKCETMVY HPNIDLEGNV CLNILREDWK PVLTINSIY GLQYLFLEPN PEDPLNKEAA EVLQNNRRLF EQNVQRSMRG GYIGSTYFER CLK
Formulation	A 0.2 µm filtered concentrated solution in 50 mM HEPES, pH 6.5, with 125 mM NaCl, 10 % Glycerol, 5 % Trehalose, 1 mM DTT.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 6 months from date of receipt, -20 to -70 °C as supplied.- 3 months, -20 to -70 °C under sterile conditions after opening.

Note: This product is for in vitro research use only