

Recombinant Human Retinoblastoma-Associated Protein Fragment, 137a.a. with His-tag



Catalog No: #AP60452

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP60452-1 10ug #AP60452-2 100ug #AP60452-3 500ug

Description

Product Name	Recombinant Human Retinoblastoma-Associated Protein Fragment, 137a.a. with His-tag
Host Species	Escherichia coli.
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Other Names	p105-Rb, pRb, pp110
Calculated MW	Approximately 16.5 kDa, a single non-glycosylated polypeptide chain containing 146 amino acids.
Target Sequence	MASFPSSPLR IPGGNIYISP LKSPYKISEG LPTPTKMTPR SRILVSIGES FGTSEKFQKI NQMVCNSDRV LKRSAEGSNP PKPLKCLRFD IEGSDEADGS KHLPGESKFQ QKLAEMTSTR TRMQKQKMND SMDTSNKEEK HHHHHH
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Background

Rb encoded by the RB1 gene in humans, is expressed by retina and belongs to the retinoblastoma-associated protein family. The whole protein consists of 928 a.a. and the rHuRb fragment occupies sequence of 792-929 a.a.. Rb is a key regulator of entry into cell division that acts as a tumor suppressor. It has many functions, for example, promotes G0-G1 transition when phosphorylated by CDK3/cyclin-C, and acts as a transcription repressor of E2F1 target genes and so on. The rHuRb is the region that rich of modified residue like phosphothreonine and N6-acetylysine.

Note: This product is for in vitro research use only