

## Recombinant Human Fatty-acid-binding Protein 2

Catalog No: #AP60474


Package Size: #AP60474-1 5ug #AP60474-2 100ug #AP60474-3 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Fatty-acid-binding Protein 2
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	P12104
GeneID	2169
Calculated MW	Approximately 15.1 kDa, a single non-glycosylated polypeptide chain containing 131 amino acids.
Target Sequence	AFDSTWKVDR SENYDKFMEK MGVNIVKRKL AAHDNLKLT I TQEGNKFTVK ESSAFRNIEV VFELGVTFNY NLADGTELRG TWSLEGNKLI GFKFRDNGN ELNTVREIIG DELVQTYVYE GVEAKRIFKK D
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

## Images



Note: This product is for in vitro research use only