

Recombinant Human Uteroglobin

Catalog No: #AP60479


Package Size: #AP60479-1 10ug #AP60479-2 100ug #AP60479-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Uteroglobin
Host Species	Escherichia coli.
Purification	> 98 % by SDS-PAGE and HPLC analyses.
Uniprot	P11684
GeneID	7356
Calculated MW	Approximately 15.8 kDa, a disulfide-linked homodimeric protein containing two 70 amino acid residues polypeptide.
Target Sequence	EICPSFQRVI ETLMDTPSS YEAAMELFSP DQDMREAGQ LKKLVDLTPQ KPRESIIKLM EKIAQSSLCN
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Images



Note: This product is for in vitro research use only