

Recombinant Rat Uteroglobin

Catalog No: #AP60498



Package Size: #AP60498-1 10ug #AP60498-2 100ug #AP60498-3 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Uteroglobin
Host Species	Escherichia coli.
Purification	> 97 % by SDS-PAGE and HPLC analyses.
Uniprot	P17559
GeneID	25575
Calculated MW	Approximately 17.0 kDa, a homodimeric protein consisting of two 77 amino acid non-glycosylated polypeptide chains.
Target Sequence	SSDICPGFLQ VLEALLGSE SNYEALKPF NPASDLQNAG TQLKRLVDTL PQETRINIVK LTEKILTSPL CEQDLRV
Formulation	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note: This product is for in vitro research use only