

# Recombinant HBV Surface Antigen (adr)

Catalog No: #AP60502



Package Size: #AP60502-1 100ug #AP60502-2 500ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant HBV Surface Antigen (adr)
Host Species	CHO
Purification	> 95 % by HPLC analyses.
Calculated MW	Approximately 23.0-27.0 kDa on a reducing SDS-PAGE.
Formulation	Sterile Filtered solution containing 20 mM PB, 154 mM NaCl.
Storage	Do not freeze.- Refer to lot specific COA for the Use by Date when stored at 2 to 8 °C as supplied.- 3 months, 2 to 8 °C under sterile conditions after opening.

Note: This product is for in vitro research use only