

Recombinant Single-stranded DNA Binding Protein

Catalog No: #AP60509


Package Size: #AP60509-1 100ug #AP60509-2 500ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Single-stranded DNA Binding Protein
Host Species	Escherichia coli
Purification	> 95 % by SDS-PAGE and HPLC analyses.
Uniprot	Q97W73
GeneID	1453842
Calculated MW	Approximately 16.1 kDa, a single non-glycosylated polypeptide chain containing 148 amino acids.
Target Sequence	MEEKVGNLKP NMESVNVTVR VLEASERQI QTKNGVRTIS EAIVGDETGR VKLTLWGKHA GSIKEGQVVK IENAWTTAFK GQVQLNAGSK TKIAEASEDG FPESQIPEN TPTAPQQMRG GGRGFRGGGR RYGRRGRRQ ENEEGEEE
Formulation	0.2 µm filtered concentrated solution in 20 mM Tris, pH 7.4, 200 mM NaCl, 50 % Glycerol, 1 mM EDTA, 0.5 mM DTT.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 6 months from date of receipt, -20 to -70 °C as supplied.- 3 months, -20 to -70 °C under sterile conditions after opening.

Images



Note: This product is for in vitro research use only