

# Recombinant Mouse Gastrin-releasing peptide receptor (Grpr), partial

Catalog No: #AP89671

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP89671-1 20ug #AP89671-2 100ug #AP89671-3 1mg

## Description

Product Name	Recombinant Mouse Gastrin-releasing peptide receptor (Grpr), partial
Brief Description	Recombinant Proteins
Host Species	E.coli
Other Names	Grpr; Gastrin-releasing peptide receptor; GRP-R; GRP-preferring bombesin receptor
Uniprot	P21729
Target Species	Mus musculus (Mouse)
Tag Info	N-terminal His-tagged Tag-Free
Target Sequence	Protein Length:Partial Expression Region: 327-384aa SKSFRKQFNTQLLCCQPGLMNRSHSTGRSTTCMTSFKSTNPSATFSLINRNICHEGYV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

&gt;85% (SDS-PAGE)

Note: This product is for in vitro research use only