

Recombinant Human Amphiphysin (AMPH)

Catalog No: #AP89672



Package Size: #AP89672-1 20ug #AP89672-2 100ug #AP89672-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Amphiphysin (AMPH)
Brief Description	Recombinant Proteins
Host Species	E.coli
Other Names	Amp; AMPH 1; AMPH; AMPH_HUMAN; AMPH1; Amphiphysin (Stiff-Man syndrome with breast cancer 128kD autoantigen); Amphiphysin (Stiff-Man syndrome with breast cancer 128kDa autoantigen); Amphiphysin; Amphiphysin I; CG8604; Damp; dAmph; Dmel\CG8604; HGNC:471; MGC49429; Stiff man syndrome with breast cancer 128kD autoantigen; Stiff man syndrome with breast cancer 128kDa autoantigen
Uniprot	P49418
Target Species	Homo sapiens (Human)
Tag Info	N-terminal His-tagged Tag-Free
Target Sequence	Protein Length:Full length protein Expression Region:1-695 MADIKTGIFA KNVQKRLNRA QEKVLQKLGK ADETKDEQFE EYVQNFKRQE AEGTRLQREL RGYLAAIKGM QEASMKLTES LHEVYEPDWY GREVDKMGVE KCDVLWEDFH QKLVDSLLT LDTYLGQFPD IKNRIAKRSR KLVYDYSARH HLEALQSSKR KDESRIKAE EEFQAQKVF EEFNVDLQEE LPSLWSRRVG FYVNTFKNVS SLEAKFHKEI AVLCHKLYEV MTKLGDQHAD KAFTIQGAPS DSGPLRIAKT PSPPEEPSPL PSPTASPNHT LAPASPAPAR PRSPSQTRKG PPVPLPKVT PTKELQENI ISFFEDNFVP EISVTTPSQN EVPEVKKEET LLDLDFDPFK PEVTPAGSAG VTHSPMSQTL PWDLWTTSTD LVQPASGGSF NGFTQPQDTS LFTMQTDQSM ICNLAESEQA PPTPEKAEPP LAAVTPAVGL DLGMDTRAEE PVEEAVIIPG ADADAAVGT LVSAAEGAPGE EAEAEKATVP AGEGVSLEEA KIGTETTEGA ESAQPEAEEL EATVPQEKVI PSVVIENASN HEEEGENEIT IGAEPKETTE DAAPPGPTSE TPELATEQKP IQDPQPTPSA PAMGAADQLA SAREASQELP PGFLYKVELT HDFEAANSDE LTLQRGDVVL VVPSDSEADQ DAGWLVGVKE SDWLQYRDLA TYKGLFPENF TRRLD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only