

# Recombinant Human Dual specificity protein phosphatase 6 (DUSP6)

Catalog No: #AP89673

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP89673-1 20ug #AP89673-2 100ug #AP89673-3 1mg

## Description

Product Name	Recombinant Human Dual specificity protein phosphatase 6 (DUSP6)
Brief Description	Recombinant Proteins
Other Names	Dual specificity protein phosphatase PYST1 (Mitogen-activated protein kinase phosphatase 3) (MAP kinase phosphatase 3) (MKP-3)
Uniprot	Q16828
Target Species	Homo sapiens (Human)
SDS-PAGE MW	49.8 kDa
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Target Sequence	<p>Protein Length: Full length</p> <p>Expression Region: 1-381aa</p> <p>MIDTLRPVPPFASEMAISKTVAWLNEQLELGNERLLLMDCRPQELYESHIESAINVAIPGIMLRRRLQKGNLPVRA            LFTRGEDRDRFRTRRCGTDTVVLYDESSDWNENTGGESVLGLLLKCLKDEGCRAFYLEGGFSKFQAEFSLHC            ETNLDGSCSSSPPLPVLGLGGLRISSDSSSDIESDLDRDPNSATDSGSPLSNSQPSFPVEILPFLYLGCARD            STNLDVLEEFGIKYILNVTPNLPNLFENAGEFKYKQIPISDHWSQNLQQFFPEAISFIDEARGKNCVGLVHCLAGI            SRSVTVTVAAYLMQKLNLSMNDAYDIVKMKKSNISPNFNFMGQLLDFERTLGLSSPCDNRVPAQQLYFTTPSNQ            NVYQVDSLQST</p>
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

Greater than 90% as determined by SDS-PAGE.

Note: This product is for in vitro research use only