

T2R5 Polyclonal Antibody

Catalog No: #48849

Package Size: #48849-1 50ul #48849-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	T2R5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	T2R5 Polyclonal Antibody detects endogenous levels of T2R5 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human TAS2R5.
Target Name	TAS2R5
Other Names	TAS2R5; Taste receptor type 2 member 5; T2R5
Uniprot	Q9NYW4
GeneID	54429
Calculated MW	35kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000. ELISA: 1:5000.

Images



Western blot analysis of the lysates from HepG2 cells using TAS2R5 antibody.

Note: This product is for in vitro research use only