

Rat Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit

Catalog No: #EK20594



Package Size: #EK20594-1 48T #EK20594-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.36ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only