

Human N6-adenosine-methyltransferase 70 kDa subunit (METTL3) ELISA Kit



Catalog No: #EK20597

Orders: order@signalwayantibody.com

Package Size: #EK20597-1 48T #EK20597-2

Support: tech@signalwayantibody.com

Description

Product Name	Human N6-adenosine-methyltransferase 70 kDa subunit (METTL3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.103ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only