

Mouse Pyruvate Dehydrogenase Kinase Isozyme 1 (PDK1) ELISA Kit

Catalog No: #EK20600

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20600-1 48T #EK20600-2 96T

Description

Product Name	Mouse Pyruvate Dehydrogenase Kinase Isozyme 1 (PDK1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.055ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only