

MBP Rabbit mAb

Catalog No: #58496

Package Size: #58496-1 50ul #58496-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MBP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	All
Specificity	MBP antibody detects endogenous levels of C-terminal, internal, and N-terminal MBP-tagged proteins.
Immunogen Description	A synthesized peptide
Other Names	MalE; malJ; Maltodextrin binding protein; Maltose ABC transporter periplasmic protein; Maltose binding periplasmic protein; MMBP; Periplasmic maltose binding protein;
Accession No.	Uniprot:P0AEX9
Uniprot	P0AEX9
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000

Images



Western blot analysis of MBP expression in E.coli lysate.

Product Description

MBP-Tag antibody is part of the Tag series of antibodies, the best quality in the research.

Background

MBP-Tag antibody is part of the Tag series of antibodies, the best quality in the research.

Note: This product is for in vitro research use only