

RFP Rabbit mAb

Catalog No: #58505

Package Size: #58505-1 50ul #58505-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RFP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	All
Specificity	RFP Antibody detects endogenous levels of total RFP
Immunogen Description	A synthesized peptide derived from human RFP
Other Names	GFP like chromoprotein; red fluorescent protein; Red fluorescent protein eqFP611; RFP;
Accession No.	Uniprot:Q8ISF8
Uniprot	Q8ISF8
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:1000~1:5000 ICC/IF 1:50~1:200

Images



Western blot analysis of RFP expression in (1) 293T cell lysate; (2) 293T cell lysate transfected with RFP.



Immunofluorescent analysis of 293 cells, using RFP Antibody

Product Description

Red fluorescent protein cloned from *Entacmaea quadricolor* eqFP611 (Uniprot IDQ8ISF8) is an important tool in biotechnology and cell biology.

Background

Red fluorescent protein cloned from *Entacmaea quadricolor* eqFP611 (Uniprot IDQ8ISF8) is an important tool in biotechnology and cell biology.

Note: This product is for in vitro research use only