

Topoisomerase II alpha Rabbit mAb

Catalog No: #58521

Package Size: #58521-1 50ul #58521-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Topoisomerase II alpha Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	Topoisomerase II alpha Antibody detects endogenous levels of total Topoisomerase II alpha
Immunogen Description	A synthesized peptide derived from human Topoisomerase II alpha
Other Names	DNA topoisomerase II; DNA topoisomerase II-alpha; TOP2; TOP2A; TP2A; alpha isozyme;
Accession No.	Uniprot:P11388
Uniprot	P11388
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of Topoisomerase II alpha expression in (1) HeLa cell lysate; (2) Jurkat cell lysate.



Immunohistochemical analysis of paraffin-embedded human colon cancer, using Topoisomerase II alpha Antibody.

Product Description

DNA topoisomerase 2-alpha; DNA topoisomerase II, 170 kD; DNA topoisomerase II, alpha isozyme; TOP2; TOP2A; topoisomerase (DNA) II alpha 170kDa; TP2A.

Background

DNA topoisomerase 2-alpha; DNA topoisomerase II, 170 kD; DNA topoisomerase II, alpha isozyme; TOP2; TOP2A; topoisomerase (DNA) II alpha 170kDa; TP2A.

Note: This product is for in vitro research use only