

Plectin Rabbit mAb

Catalog No: #58595

Package Size: #58595-1 50ul #58595-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Plectin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	Plectin Antibody detects endogenous levels of total Plectin
Immunogen Description	A synthesized peptide derived from human Plectin
Other Names	EBS1; EBSO; HD1; Hemidesmosomal protein 1; PCN; pleC; PLEC1 ;Plectin 1; Plectin 6; Plectin; PLTN
Accession No.	Uniprot:Q15149
Uniprot	Q15149
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:2000 IHC 1:50~1:100 IP 1:30

Images



Western blot analysis of Plectin expression in Jurkat cell lysate.

Product Description

Interlinks intermediate filaments with microtubules and microfilaments and anchors intermediate filaments to desmosomes or hemidesmosomes. Could also bind muscle proteins such as actin to membrane complexes in muscle. May be involved not only in the filaments network, but also in the regulation of their dynamics.

Background

Interlinks intermediate filaments with microtubules and microfilaments and anchors intermediate filaments to desmosomes or hemidesmosomes. Could also bind muscle proteins such as actin to membrane complexes in muscle. May be involved not only in the filaments network, but also in the regulation of their dynamics.

Note: This product is for in vitro research use only