

Hsp90 alpha Rabbit mAb

Catalog No: #58718

Package Size: #58718-1 50ul #58718-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Hsp90 alpha Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Hsp90 alpha antibody detects endogenous levels of total Hsp90 alpha
Immunogen Description	A synthesized peptide derived from human Hsp90 alpha
Other Names	FLJ31884; HSPC1 ; EL52; Hsp90; HSPCA; LAP2 ; HSPN ; HSP90A; HSP 86; heat shock 90kD protein 1 alpha antibody;
Accession No.	Uniprot:P07900
Uniprot	P07900
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200


Images



Western blot analysis on HeLa cell lysate using HSP90A Antibody.



Immunohistochemical analysis of paraffin-embedded human breast cancer, using Hsp90 alpha antibody.



Immunofluorescent analysis of HeLa cells, using Hsp90 alpha Antibody.

Product Description

HSP90A a molecular chaperone of the heat shock protein 90 family. Has ATPase activity. Known to interact with a wide variety of proteins including steroid hormone receptors, neuropeptide Y, FKBP51/54, and FKBP52. G protein-coupled receptor kinases are stabilized by interacting with HSP 90. Hsp70 and Hsp90 promote tau solubility and tau binding to microtubules, reducing insoluble tau.

Background

HSP90A a molecular chaperone of the heat shock protein 90 family. Has ATPase activity. Known to interact with a wide variety of proteins including steroid hormone receptors, neuropeptide Y, FKBP51/54, and FKBP52. G protein-coupled receptor kinases are stabilized by interacting with HSP 90. Hsp70 and Hsp90 promote tau solubility and tau binding to microtubules, reducing insoluble tau.

Note: This product is for in vitro research use only