

YAP1 Rabbit mAb

Catalog No: #58764

Package Size: #58764-1 50ul #58764-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	YAP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	YAP1 Antibody detects endogenous levels of total YAP1
Immunogen Description	A synthesized peptide derived from human YAP1
Other Names	65 kDa Yes-associated protein; YAP1; YAP65;
Accession No.	Uniprot:P46937
Uniprot	P46937
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:5000~1:10000 IHC 1:50~1:100 ICC/IF 1:50~1:100 IP 1:50 FC 1:50


Images



Western blot analysis of YAP expression in HeLa cell lysates.



Immunohistochemical analysis of paraffin-embedded human uterus, using YAP1 Antibody.



Immunofluorescent analysis of MCF7 cells, using YAP1 Antibody .

Product Description

Plays a key role to control cell proliferation in response to cell contact. Phosphorylation of YAP1 by LATS1/2 inhibits its translocation into the nucleus to regulate cellular genes important for cell proliferation, cell death, and cell migration.

Background

Plays a key role to control cell proliferation in response to cell contact. Phosphorylation of YAP1 by LATS1/2 inhibits its translocation into the nucleus to regulate cellular genes important for cell proliferation, cell death, and cell migration.

Note: This product is for in vitro research use only