

## SOX18 Rabbit mAb

Catalog No: #58793

Package Size: #58793-1 50ul #58793-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	SOX18 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	SOX18 Antibody detects endogenous levels of total SOX18
Immunogen Description	A synthesized peptide derived from human SOX18
Other Names	HLTS; Sox18 ; SRY box 18; Transcription factor SOX 18;
Accession No.	Uniprot:P35713
Uniprot	P35713
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200

## Images



Western blot analysis of SOX18 expression in (1) K562 cell lysate; (2) human stomach lysate.

## Product Description

Binds to the consensus sequence 5'-AACAAAG-3' and is able to trans-activate transcription via this site. Defects in SOX18 are the cause of hypotrichosis-lymphedema-telangiectasia syndrome (HLTS) [MIM:607823].

## Background

---

Binds to the consensus sequence 5'-AACAAAG-3' and is able to trans-activate transcription via this site. Defects in SOX18 are the cause of hypotrichosis-lymphedema-telangiectasia syndrome (HLTS) [MIM:607823].

---

Note: This product is for in vitro research use only