

NUP98 Rabbit mAb

Catalog No: #58817

Package Size: #58817-1 50ul #58817-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NUP98 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	NUP98 Antibody detects endogenous levels of total NUP98
Immunogen Description	A synthesized peptide derived from human NUP98
Other Names	ADAR2; Nup98; Nucleoporin Nup98; NUP96; Nucleoporin Nup96; ADIR2;
Accession No.	Uniprot:P52948
Uniprot	P52948
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of NUP98 expression in HeLa cell lysate.



Immunohistochemical analysis of paraffin-embedded human testis, using NUP98 Antibody.

Product Description

Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. Nup98 and Nup96 are involved in the bidirectional transport across the NPC. May anchor NUP153 and TPR to the NPC.

Background

Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. Nup98 and Nup96 are involved in the bidirectional transport across the NPC. May anchor NUP153 and TPR to the NPC.

Note: This product is for in vitro research use only